

MAGOHB-mago-nashi homolog B (Drosophila) Gene

Reference: PTXBC010905

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MGN2; mago; magoh; protein mago nashi homolog 2; mago-nashi homolog B; mago homolog B, exon junction complex core component		

For research use only.

GENE DESCRIPTION	mago-nashi homolog B (Drosophila)
GENE ID	55110
ORIGIN SPECIES	Human
ACCESSIONS	BC010905
VECTOR	pDONR223
NCBI SYMBOL	MAGOHB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtggctagcgatttctacctgcgctactacgttagggcacaaggcgaagtgtgggcacagtttctggagttcgaattcggccgacggaaagcttagatgccaacaacagcaattacaaaaatgatgtgatgacagaaagaggcttatgtgcacaagagtgaatgaaga
actgaagagaattattgatgacagtgaaattacaaaagaagatgatgctttgtggcctccccctgataggggtggccgacaggagctgaaattgtaattggagatgagcacatatctttaccacatcaaaaataggttctctattgatgtaaatcagtcgaaagatcctgaaggccttcgagtat
tttactatttgggtacaagactgaaatgttttagtttcagcttattggattacacttcaagattaaccaatttaa