

HBG1-hemoglobin, gamma A Gene

Reference: PTXBC010914

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HBG-T2; HBGA; HBGR; HSGGL1; PRO2979; hemoglobin subunit gamma-1; A-gamma globin; gamma A hemoglobin; gamma globin; gamma-1-globin; hb F Agamma; hemoglobin gamma-1 chain; hemoglobin gamma-a chain; hemoglobin, gamma A; hemoglobin, gamma, regulator of; hemog		

For research use only.

GENE DESCRIPTION	hemoglobin, gamma A
GENE ID	3047
ORIGIN SPECIES	Human
ACCESSIONS	BC010914
VECTOR	pDONR223
NCBI SYMBOL	HBG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtcattcacagaggagacaaggctactatcacaagcctgtggggcaaggtgaatgtggaagatctggaggagaacacctgggaaggtcctgtgttctacccatggaccagaggtctttgacagctttggcaacctgtcctctgacctcatgggcaacccaaagt
caaggcacatggcaagaaggtgctgacttcttggagatgccataaagcacctggatgatctcaaggcaccctttgccagctgagtgaaatgcaatgacatgtggaatctcaagctcctgggaaatgtgctgtgacctgtttggcaatcatttcggcaaga
attcacccctgaggtgcaggtctcctggcagaagatggtgactggagtgccagtgccctgtctccagataccactga