

RPL35-ribosomal protein L35 Gene

Reference: PTXBC010919

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L35; 60S ribosomal protein L35; ribosomal protein L35		
For research use only.			

GENE DESCRIPTION	ribosomal protein L35
GENE ID	11224
ORIGIN SPECIES	Human
ACCESSIONS	BC010919
VECTOR	pDONR223
NCBI SYMBOL	RPL35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaagatcaaggctcagatctcgcgggaagaagaaggaggagctgaaacagctggacgacctgaagggtggagctgtccagctgcgcgtcgccaaagtgcacaggcgggtgcggcctccaagctcttaagatccgagtcgtccggaatccattgccgtgttctcacag
ttattaccagactcagaagaaaacctcaggaaattctacaagggaagaagtacaagccctggacgtgcgcctaagaagacacgtgcatgcgcggcggctcaacaagcagcaggagaacctgaagaccaagaagcagcagcgggaaggagcggctgtacctgctgcggaa
gtacgcgggtcagggcctga