

CLDN10-claudin 10 Gene

Reference: PTXBC010920

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPETRL3; OSP-L; claudin-10; OSP-like protein; claudin 10		
For research use only.			

GENE DESCRIPTION	claudin 10
GENE ID	9071
ORIGIN SPECIES	Human
ACCESSIONS	BC010920
VECTOR	pDONR223
NCBI SYMBOL	CLDN10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctagcacggcttcggagatcatgccttcattgggtccatctcaggctgggtggtgtctccacgtgccaccgactactggaaggtgtctaccatcgacggcagcgtcatcacaccgccacctattggccaacctgtggaaggcgtgcgttacgactccacggcgctc
caactgcaaggacttccctccatgctggcgctggacggttatatacaggcatgtagaggacttatgatcgtctgcagcctgggtctcttgggtccatattgcgctcttggatgaagtgtaccaaagtcggaggtccgataaaagccaaagctaaattgctgttggctggatgtat
tcatactgtcagggtgtgtcgaatgactggatgtccctatatgcaacaaaaacacacgaattctttgatcctctcttggagcaaaagtatgaattaggagccgctctgtttatggatggcaggaagccctcactgtgcataattgggtgtgtcatatttgcctttcaatatctgacaacaac
aaaacaccagatacacatacaacggggccacatctgtcatgtcttcggacaaagtatcgtgggagaagattttaaacaacaaccctcaaacagtttgataaaatgcttatgtctaa