

HLA3-HERV-H LTR-associating 3 Gene

Reference: PTXBC010922

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HERV-H LTR-associating protein 3; human endogenous retrovirus-H long terminal repeat-associating 3; human endogenous retrovirus-H long terminal repeat-associating protein 3; HERV-H LTR-associating 3		

For research use only.

GENE DESCRIPTION	HERV-H LTR-associating 3
GENE ID	11147
ORIGIN SPECIES	Human
ACCESSIONS	BC010922
VECTOR	pDONR223
NCBI SYMBOL	HLA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcggtgcatgttataaacaaccactaaacctccggatctgagcctcccgagagggaatgcagaatgacgccagggcagcagatgtgatgtcaccgagatgcagaggatactcagtcaccaacattactgagcatctactctgtgccgtatgtctgtcaacggaaaggggtcc
ctatccagacccaagagagcattcttggatctcttgcagaagaatttgaggcgaatccatagagtaagcttagtgatgtgtcagacctctgagcccaagcaaaagccatcatatccctgtgacctgcatgtatataccagatggcctgaagcaagtgaagaatcacaagaagtga