

LDB3-LIM domain binding 3 Gene

Reference: PTXBC010929

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMD1C; CMH24; CMPD3; CYPHER; LDB3Z1; LDB3Z4; LVNC3; MFM4; ORACLE; PDLIM6; ZASP; LIM domain-binding protein 3; PDZ and LIM domain 6; Z-band alternatively spliced PDZ-motif protein; cardiomyopathy, dilated 1C (autosomal dominant); protein cypher; LIM domain		

For research use only.

GENE DESCRIPTION	LIM domain binding 3
GENE ID	11155
ORIGIN SPECIES	Human
ACCESSIONS	BC010929
VECTOR	pDONR223
NCBI SYMBOL	LDB3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttacagtgtgacctgactgggcccgggcccgtgggcttcgtctgaggggggcaaggacttcaacatgccctcactatctcccggatcacaccaggcagcaaggcagccagtcagctcagccaggggtgacctgtgtgtggccattgacggcgtcaacacagacacca
tgaccacctggaagcccaagaacaagatcaagtcgagctacaaactgagcctcaccctgcagaaatcaaaagctcccaattcccatctccacgacagcactccagtcagaccctctgcccgtgatccctcaccagaagggtgtagtcaactctccagccaacgccgactaccagg
aacgttcaacccagtgccctgaaggactcggccctgtccaccacaagcccatcgaggtgaaggggctggcgccgaaggcccatcatccatgcgcagtcacaacagcccatcagcatgtattccaggaatgcatcatgagccatcgctgggcagggccaagccaaggca
gtgacttcagtgaggcctccctattaaggaccttgcgtagacagcgctccccgtctaccaggctgtgattaagagccagaacaagccagaagatgaggtgacgagtggtggcagccgttctcctcaacctgcagtcctcctccgcatctggccagatgacgggacagaat
tcacgcaagaccctgatgaagaagctctgcgaaggtcaagggaaaggttgaacggaacgtaacagccacggtttgccaattgcgcaactggcaccatggccttcagcccaaatccttaattgtaaaagctaa