

RWDD2A-RWD domain containing 2A Gene

Reference: PTXBC010930

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RWDD2; dJ747H23.2; RWD domain-containing protein 2A; RWD domain containing 2; RWD domain containing 2A		

For research use only.

GENE DESCRIPTION	RWD domain containing 2A
GENE ID	112611
ORIGIN SPECIES	Human
ACCESSIONS	BC010930
VECTOR	pDONR223
NCBI SYMBOL	RWDD2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgcttgcgtgaagaaagccttcagcttgcagctactggagatggaaatgctgtttctatgtttcctaaccagagagaagtaaaacttgaagatgtaaatgccctgacgaatataagaggtatttgaaggcacaaggaggcgctgccacaaaaatcgaattgtaattacactcca
gattgaagagcccaagggtgaaaattgatttgcagtgaccatgcctcacagctacccctatttagcattgcagctgttttgacggtcatctgaacttgacagacatcagcagctacttctcaacaaggctcacttctacatagggaactttgatccagggtgagctctgtgtatgtgcagcaatcc
agtggctacaggacacagctgcatcttatttctgaacagaagcttgcataatgaacatctacacaagcaagccagctcaagaacacatctccgaatgtggtatctacagtcaccatataatcagcaggaacaaaggaataatttggatgttgggaaaaggttagatgtgactggatttt
gcatgacaggaagcctggtataatctgtgtggaggtttcaaaagcactgtgaggagttctggcacacaattaggtaccccaattggaagcacatttctgtaagcatgctgaaagcgtggagacagaaggaatggtgaggacctgcgcttttcattctttgaagagttactcctga
ggctcatggtgactatggatlaaggaatgactatcacatgaatctggcccaattctagaatttctcaagaagcacaagaagtgcagcatgttttcagatactatttggatgaaagcaaaagttcagactcataa