

CSN3-casein kappa Gene

Reference: PTXBC010935

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CNS10; CSN10; CSNK; KCA; kappa-casein; casein kappa		
For research use only.			

GENE DESCRIPTION	casein kappa
GENE ID	1448
ORIGIN SPECIES	Human
ACCESSIONS	BC010935
VECTOR	pDONR223
NCBI SYMBOL	CSN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaaagtgtttcttagtgcacatgcctggcattaacctgctcttttgggtgtggaggttcaaaaccagaaacaccagcatgccaatgagaatgatgaagaccattctatcagaaaacagctccatgtgccaatgtattatgccaatagctatccttattatgaaccaattgtacc
aacgtagaccagctatagcaataataatccatgtgtcctgcacatatattatgcaaacccagctgtagttaggccacatgcccaatctctcagcggcaataacctgccaaatagccaccaccctgtggtacgtgcccaaacctgcatccatcatttattgcatcccccaagaaaat
caggataaaaataatcatccctaccatcaataccattgtactgttgaacctacaccagctcctgccactgaaccaacgggtggacagtgtagtcactccagaagctttttcagagtgccatcatcacgagcaccctgagacaaccacagttgcagttactccacacacgacataa