

SRP19-signal recognition particle 19kDa Gene

Reference: PTXBC010947

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	signal recognition particle 19 kDa protein; signal recognition particle 19kD; signal recognition particle 19kDa; signal recognition particle 19		

For research use only.

GENE DESCRIPTION	signal recognition particle 19kDa
GENE ID	6728
ORIGIN SPECIES	Human
ACCESSIONS	BC010947
VECTOR	pDONR223
NCBI SYMBOL	SRP19
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgcactccgcgcgtcccccggccgaccaggacaggtttattgtatctatcctctattttaataataagaaccatcgagaggggaaggcgaatcccataagtaaggctgttgaataatctacagctacagagattcaagatgtatgttcagcagttggacttaacgtattctt
gagaaaaataaatgtactctagagaatggaatcgtgatgtccaatacagaggcagagtcgggtccagctcaaacaggaaagatgggagcctctgcttgcagttccatcacgtaagtcagtaattgtgtatgcagcagaaatgatactaaactaaaacaaggacacaaaaacag
gaggtgctgaccaaagtcttcaacaaggaggggaagtaaaaaagggaaggaaagaaaaagaagtaa