

PIP-prolactin-induced protein Gene

Reference: PTXBC010951

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GCDFP-15; GCDFP15; GPIP4; prolactin-inducible protein; SABP; gross cystic disease fluid protein 15; secretory actin-binding protein; prolactin induced protein		

For research use only.

GENE DESCRIPTION	prolactin-induced protein
GENE ID	5304
ORIGIN SPECIES	Human
ACCESSIONS	BC010951
VECTOR	pDONR223
NCBI SYMBOL	PIP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgcttgcctcagctcctgttcaggccagccctgcctcctgttctctgcctgcagttgggggccaacaaagtcaggacacactcggagatcataataagaattttgacattccaagtcagtcacaaatgacgaagtcagtcagtcgttcagttcaaacagaat
tgaagaatgcatgggtggttaaaacttacctcattagcagcatccctctacaaggctgcatftaactataagtatactgcctgcctatgtgacgacaatccaaaaccttctactgggacttttacccaacagaactgtgcaaatgcagccgtcgttgatgtattcggaattagcctcgcct
gatgatgctgctgtaatccccatcaaaaacaatcggttttatactattgaaatcctaaaggtagaataa