

FGF7-fibroblast growth factor 7 (keratinocyte growth factor) Gene

Reference: PTXBC010956

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HBGF-7; KGF; fibroblast growth factor 7; FGF-7; heparin-binding growth factor 7; keratinocyte growth factor		

For research use only.

GENE DESCRIPTION	fibroblast growth factor 7 (keratinocyte growth factor)
GENE ID	2252
ORIGIN SPECIES	Human
ACCESSIONS	BC010956
VECTOR	pDONR223
NCBI SYMBOL	FGF7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacaaatggatactgacatggatcctgccaactttgctctacagatcatgctttacattatctgtctagtgggtactatcttttagcttgaatgacatgactccagagcaaatggctacaaatgtgaactgttccagccctgagcgcacacacaagaagttatgattacatggaaggagagg
atataagagtgagaagactcttctgtcgaacacagtggtacctgaggatcgataaaagaggcaagtaaaaggacccaagagatgaagaataattacagtaagtaa