

CCND2-cyclin D2 Gene

Reference: PTXBC010958

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAK0002; MPPH3; G1/S-specific cyclin-D2		
For research use only.			

GENE DESCRIPTION	cyclin D2
GENE ID	894
ORIGIN SPECIES	Human
ACCESSIONS	BC010958
VECTOR	pDONR223
NCBI SYMBOL	CCND2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctgctgtgccacgaggtggaccgggtccgacggggcgtgctgggacccgaacctgctccgagacgacgctcctgcagaacctgctaccatcgaggagcgtaccttcgcagtgctcctactcaagtgcgtgcagaaggacatccaacctacatgcgcagaatggtg
gccacctggatgctggaggtctgtgaggaacagaagtcggaagaaggtcttccctctggccatgaattacctggaccgtttctggctgggggtccgactccgaagtcctatcgcaactcctgggtgctgtgcatgttctggctccaaactcaaaagaccagcccgtgaccgc
ggagaagctgtgcatcaccgacaactccatcaagccctcaggagctgctggagtggaactggtgtgctggggaagttgaagtggaacctggcagctgtcactcctcatgactcattgagcacatcttgcgaagctgccccagcagcgggagaagctgtctctgacccgaagcat
gctcagacctcattgctgtgtgcccacggactttaagtttgccatgtaccacccgtcgtatgcgaactggaagtggtggagcagccatctgtgggtccagcaggatgaggaagtgagctcgtcactgtgtatgccctgactgagctgtgggtaagatcaccacacagacgtggat
tgtctcaagcttgcaggagcagattgaggcgggtgctcctcaatagcctgcagcagtcaccgcaggaccaacgtgacgggatccaagtcggaggtgaactggaccaagccagcaccctacagacgtgcgggatatgcacctgtga