

BNIP1-BCL2/adenovirus E1B 19kDa interacting protein 1 Gene

Reference: PTXBC010959

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NIP1; SEC20; TRG-8; vesicle transport protein SEC20; BCL2/adenovirus E1B 19 kDa protein-interacting protein 1; BCL2/adenovirus E1B 19kDa interacting protein 1; transformation-related gene 8 protein; BCL2 interacting protein 1		

For research use only.

GENE DESCRIPTION	BCL2/adenovirus E1B 19kDa interacting protein 1
GENE ID	662
ORIGIN SPECIES	Human
ACCESSIONS	BC010959
VECTOR	pDONR223
NCBI SYMBOL	BNIP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggctcccaagacgtccacgtccggtatgttaaccaagagattgtcaaaattgacctggaggtgaaggcgcttattcaggatattcgtgattgttcaggacccttaagtgtcttactgaactgaataactaaagataaaagagaattcaacagttgcgtcacagaatacaggacctgga
gcagttggctaaagagcaagacaagaatcagagaacaacttctactccagggaagtggagaatcacaaaaagcagatgtcagcaatcaggcctcatggaggaaagctaatctcactgcaaaatfgcaatcgacaatctagagaagcagaacttctcaggaggagatctcttaag
gcaaaggaaaaccaccaagagagcctggccagacatccagtaccatcactgagagcctcatgggatcagcaggatgatggccagcaggtccagcagagcgaggagccatgcagtctctagtcacttcttcacgaacgatcctggatgcaaatgaagaatttaagtcctatg
ggcaccatccagctggccggagcttatcacaataacaatcgccgggagctgacggacaagcttctcattcttgcgtagccctgtttctgtacggctcttatattgtgaaaaagcggctcttccattttgtga