



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TNNT1-troponin T type 1 (skeletal, slow) Gene

Reference: PTXBC010963

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ANM; NEM5; STNT; TNT; TNTS; troponin T, slow skeletal muscle; slow skeletal muscle troponin T; troponin T type 1 (skeletal, slow); troponin-T1, skeletal, slow; troponin T1, slow skeletal type		

For research use only.

GENE DESCRIPTION	troponin T type 1 (skeletal, slow)
GENE ID	7138
ORIGIN SPECIES	Human
ACCESSIONS	BC010963
VECTOR	pDONR223
NCBI SYMBOL	TNNT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctggacaccgaggagcaggaataatgagaggagcagccggaaaggaggctctggaggaggaggaggaagcccccgaagagcggcggcgggtggcagagccagaagaggaacgcccccacaaacagccgccctgtgtctctcttgcctccgcaaaagatcccaaa
ggggagcgcgttgactctgatgacatccaccgcaagcgcgatggagaaagacctgctggagctgcagacactcatcgatgtacatttcgagcagcgggaagaaggaggagagagcgtgtgtgctctgaaggagcgcattgagcggcgccggctcagagagagccggaacacgcgtt
cagaactgagaaggaacccgaacgtcagctcgaagtcaggcggagagagaagatgagaaggaagagaggaagagcccaagaagcggcgagaggatgatgccaaagaaaagaaggtctgtccaacatggggcccaatttggcgcctactgtgtcaagcagaaacaagaagcgtgtgta
agcggcagcagggcggggagatgaaaggtgtgcgactcctctccgacgctgaagaagcctctgcagcaattgactactggggcaggaagcctccggcgaagaagcccgagctctggagcttgagaagttgcacgtgcgcgcgaagctgaacagcagaa
atatgatatcaacctgtctgtacaaccctatcagcccccagaagatttccgcgaagctggcagagcgaagccgcctctgtgagccgctggaagaatga

More info: Please contact us for any question: contact@proteogenix.fr