



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CAB39L-calcium binding protein 39-like Gene

Reference: PTXBC010993

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MO25-BETA; MO2L; bA103J18.3; calcium-binding protein 39-like; MO25beta; U937-associated antigen; antigen MLAA-34; mo25-like protein; sarcoma antigen NY-SAR-79; calcium binding protein 39 like		

For research use only.

GENE DESCRIPTION	calcium binding protein 39-like
GENE ID	81617
ORIGIN SPECIES	Human
ACCESSIONS	BC010993
VECTOR	pDONR223
NCBI SYMBOL	CAB39L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr