



Reference: PTXBC010996

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HZF1; KOX9; zinc finger protein 16; zinc finger protein KOX9; zinc finger protein Kruppel type 9		

For research use only.

GENE DESCRIPTION	zinc finger protein 16
GENE ID	7564
ORIGIN SPECIES	Human
ACCESSIONS	BC010996
VECTOR	pDONR223
NCBI SYMBOL	ZNF16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]