

GLB1L3-galactosidase, beta 1-like 3 Gene

Reference: PTXBC011001

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	beta-galactosidase-1-like protein 3; galactosidase beta 1 like 3		
For research use only.			

GENE DESCRIPTION	galactosidase, beta 1-like 3
GENE ID	112937
ORIGIN SPECIES	Human
ACCESSIONS	BC011001
VECTOR	pDONR223
NCBI SYMBOL	GLB1L3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagtcctccgccccctcagccggtgtctctctggaagagaatggcgggcatcttttctgccatttatctcaggttttctcctcggttaagcagggaagagaacttcattgcttggaaagagcgcacccgctccagcctaggtttaattggtctcttgacccctctggagctgaaga
atcgatctgtgggacttggactgaaagcacaggtcgggtaagccccacttcacactggaggccacaagtctctggttcggttcctccactatttccgggtgccaggagtagtgaggagccgctgtgagctgaagcctgtggctcaatactgcaccacctatgtt
ccgtggaacctgcatgagccagaagagcgaatttgactctctgggaacctggacctggagccttcgtctgatggccgagagatcggtgtgtgtgtgattctgctccaggccctacatctgcatgagatgacctcggggcttggccagctggctcctgcaagaccccg
gttactgttgaggacaaccaagagcttcattgaagcagttgagaagtatttgaccacctgattccagagtgattcctctccagtaaccgagggcctgtcatcgcggtgcaagtggagaatgagtaggttcattcaataagataaacatacatgccgtatctccacaaggcc
ctgctgagaagagggttggagcttctctgacctctgatggtgagaacatgtgctgagtgccacacaaaggagtggtggcccatcaattgcaaaaactcaccaggatacttcaatcagcttcataaagtcaggagccagtagatgcaggactgctctcagatgcctcac
aagcagcaacttctctcttatag