

TMEM216-transmembrane protein 216 Gene

Reference: PTXBC011010

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC244; transmembrane protein 216		
For research use only.			

GENE DESCRIPTION	transmembrane protein 216
GENE ID	51259
ORIGIN SPECIES	Human
ACCESSIONS	BC011010
VECTOR	pDONR223
NCBI SYMBOL	TMEM216
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctctctttatcttggaattgaagtaattgcctgttttttggtacaaagggaacctctgccagcgaagatgccactcagttattagcgtggccttgaccttcccatctgccatgatggcctctattacctgctgctgcagacctacgtactccgctggaagccatcatgaatggcatcttg
ctcttctctgtggtcagagcttttacttgagggtcctcacttggctgcttttccagtatggacacgattga