



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CLEC3B-C-type lectin domain family 3, member B Gene

Reference: PTXBC011024

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TNA; tetranectin; plasminogen kringle 4-binding protein; tetranectin (plasminogen-binding protein); C-type lectin domain family 3 member B		

For research use only.

GENE DESCRIPTION	C-type lectin domain family 3, member B
GENE ID	7123
ORIGIN SPECIES	Human
ACCESSIONS	BC011024
VECTOR	pDONR223
NCBI SYMBOL	CLEC3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]