

SRI-sorcin Gene

Reference: PTXBC011025

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CP-22; CP22; SCN; V19; 22 kDa protein; H_RG167B05.1; calcium binding protein amplified in mutlidrug-resistant cells		

For research use only.

GENE DESCRIPTION	sorcin
GENE ID	6717
ORIGIN SPECIES	Human
ACCESSIONS	BC011025
VECTOR	pDONR223
NCBI SYMBOL	SRI
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtaccggggcctcctggcggcggcggggtactaccaggcgggtatggagggtcccgaggggcctgctttccggacaaactcaggatccgctgtatggttactttgctgctgtagctggacaggatgggcagatagatgctgatgaattgcagagatgtctgacacagtcgtgcatgctggaggatacaaaccttttaacctggagactgcccgttattggttcaatgctgtagatagatagatgctgacacaaatgggtttcaatgaattaaagaactctgggctgactgaatggctggagacaacactttatcagttttgacactgacaggagtggaacagtagaccacaagaattgcagaaggccctgacacaaatgggatttaggttagtccccaggctgtgaattcaattgcaaaacgatacagcaccaatggaaagatcaccttcgacgactacatcgccctgctgcgtcaactgagggtctttacagacagcttcgaagacgggatactgctcagcaagggtgtgtgaattcccatatgatgattcattcaatgtgcatgaggtttaa