



Reference: PTXBC011171

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1IN; C1INH; C1NH; HAE1; HAE2; plasma protease C1 inhibitor; C1 esterase inhibitor; C1-inhibiting factor; complement component 1 inhibitor; serine/cysteine proteinase inhibitor clade G member 1; serpin G1; serpin peptidase inhibitor clade G member 1; serp		

GENE DESCRIPTION	serpin peptidase inhibitor, clade G (C1 inhibitor), member 1
GENE ID	710
ORIGIN SPECIES	Human
ACCESSIONS	BC011171
VECTOR	pDONR223
NCBI SYMBOL	SERPING1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]