

SMNDC1-survival motor neuron domain containing 1 Gene

Reference: PTXBC011234

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SMNR; SPF30; TDRD16C; survival of motor neuron-related-splicing factor 30; 30 kDa splicing factor SMNrp; SMN-related protein; splicing factor 30, survival of motor neuron-related; survival motor neuron domain-containing protein 1; tudor domain containing		

For research use only.

GENE DESCRIPTION	survival motor neuron domain containing 1
GENE ID	10285
ORIGIN SPECIES	Human
ACCESSIONS	BC011234
VECTOR	pDONR223
NCBI SYMBOL	SMNDC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcagaggatttagcaaacgagctggcaagctacaaagctcagctccagcaagttgaagctgattatctggaatggagaaatgaagattgctaaaattgaagaagattacaagaattatagaactaaccaaaccttctgtcaactcaaccttctgagacgcttgcaggttcag
acagttttgcttctactcaacctactcattcatggaaagtaggagacaagtgatgagctggtgagtggaatggacagtggtatgaagcggagattgaggagatagatgaagaaatggcaccgctgcaatcacctttgctggttatggcaatgctgaagtactccactgttgaacctcaag
cctgtagaagaaggagaaaggcaaggagagacagtgccacaacccatgtcaaaaaaagaagaatgattgccagcagcgtgaatataaaaaagaagaagctttgaaaaagctcagagaataaaaagaacttgagcaggaaagagaggaccagaagtgaaatggcaacaattcaa
caacagagcctattctaaaaacaaaaagccaggtaaagaggagtttttcttcacctgagagtgtagctgttaaagttggagtaggaacctgtggaattgctgataaacctatgacacaatatcaagatacctctaataacaatgcaggcatttgatgcctcaataa