

PSD3-pleckstrin and Sec7 domain containing 3 Gene

Reference: PTXBC011238

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EFA6D; HCA67; PH and SEC7 domain-containing protein 3; ADP-ribosylation factor guanine nucleotide factor 6; epididymis tissue protein Li 20mP; exchange factor for ADP-ribosylation factor guanine nucleotide factor 6; hepatocellular carcinoma-associated ant		

For research use only.

GENE DESCRIPTION	pleckstrin and Sec7 domain containing 3
GENE ID	23362
ORIGIN SPECIES	Human
ACCESSIONS	BC011238
VECTOR	pDONR223
NCBI SYMBOL	PSD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttggttaatggagagaattttgagttcacttaataatttcccatcggtcgccataataaagcttcaggcgctcctagaagagtcgcccaaggctcgattaaggaccacactgcaggtctgaggtcactgctctgagtcctgaacaccagagccctgcagagagtggtgataaacac
atcatctctgcaagagggaacctctccccggccgcacttcactcaggcttctactgagcagcaaggacagcctgggtttcaaatgccacttccctgctttagggtccaggtgtctctgatagcgtga