



Reference: PTXBC011355

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TFIIIC63; TFIIICepsilon; TFIIIC2-63; general transcription factor 3C polypeptide 5; TF3C-epsilon; TFIIIC 63 kDa subunit; general transcription factor IIIC, polypeptide 5, 63kDa; transcription factor IIIC 63 kDa subunit; transcription factor IIIC subunit e		

For research use only.

GENE DESCRIPTION	general transcription factor IIIC, polypeptide 5, 63kDa
GENE ID	9328
ORIGIN SPECIES	Human
ACCESSIONS	BC011355
VECTOR	pDONR223
NCBI SYMBOL	GTF3C5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr