

TRAPPC2L-trafficcking protein particle complex 2-like Gene

Reference: PTXBC011369

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC176; trafficking protein particle complex subunit 2-like protein; hematopoietic stem/progenitor cells 176; trafficking protein particle complex 2 like		

For research use only.

GENE DESCRIPTION	trafficking protein particle complex 2-like
GENE ID	51693
ORIGIN SPECIES	Human
ACCESSIONS	BC011369
VECTOR	pDONR223
NCBI SYMBOL	TRAPPC2L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggtgcatcgcggtgattgccaaggagaattacccctctacattcgagcaccctacggagacgagctgaagttccactacatggtgcacacatctctggacgtggtggatgagaagatccgcaatggggaaggccctggtcgaccagaggagctgtacctgggcct
gctctacccacggaggactacaaggtatacggctacgtcaccaactccaaggtgaagtttgcattggtgtagattcccaacacagccctcgagacaacgaattcgagcatgttcgggaagctacacaactctacacagacgtgatgtgcaacccctctacaacccgggggac
cgcatccagtcaggccttgataacatggtgacgtcgatgatgatacaggtgtgctga