



Reference: PTXBC011371

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACHRB; CHRNb; CMS1D; CMS2A; CMS2C; SCCMS; acetylcholine receptor subunit beta; acetylcholine receptor, nicotinic, beta 1 (muscle); cholinergic receptor, nicotinic beta 1; cholinergic receptor, nicotinic, beta 1 (muscle); cholinergic receptor, nicotinic, b		

For research use only.

GENE DESCRIPTION	cholinergic receptor, nicotinic, beta 1 (muscle)
GENE ID	1140
ORIGIN SPECIES	Human
ACCESSIONS	BC011371
VECTOR	pDONR223
NCBI SYMBOL	CHRNA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr