

CFL2-cofilin 2 (muscle) Gene

Reference: PTXBC011444

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NEM7; cofilin-2; cofilin 2 (muscle); cofilin, muscle isoform; nemaline myopathy type 7; cofilin 2		

For research use only.

GENE DESCRIPTION	cofilin 2 (muscle)
GENE ID	1073
ORIGIN SPECIES	Human
ACCESSIONS	BC011444
VECTOR	pDONR223
NCBI SYMBOL	CFL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctctggagttacagtgaatgatgaagtcataaagttttatgatagaaagtaaggaaatctctacacagaggagatcaaaaagagaagaaagcagttctctctgttaagcgtgacaaaagacaataattgtagaggaagcaagcagatctgggtgggtgacattggtga
tactgtagaggaccctacacatctttgtgaagttgctacctctgaatgattgccgatagctttgtacgatgccacatacgaacaaaagagcttaagaaaagacacagtagtatttatactctggctctgaaagtcaccttfaaaagcaagatgattatgctagctctaaagatgccattaa
aaagaaattacaggtattaaacatgagtggaagtaaatggcttgatgataaaggaccgttcgacacttggagagaaattgggaggcaatgtagtagttcacttgaaggaaaaccattataa