

RPL32-ribosomal protein L32 Gene

Reference: PTXBC011514

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L32; PP9932; 60S ribosomal protein L32; ribosomal protein L32		

For research use only.

GENE DESCRIPTION	ribosomal protein L32
GENE ID	6161
ORIGIN SPECIES	Human
ACCESSIONS	BC011514
VECTOR	pDONR223
NCBI SYMBOL	RPL32
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgcccctcagacccttgtgaagcccaagatcgtcaaaaagagaaccaagaagttcatccggcaccagtcagaccgatatgtcaaaattaagcgttaactggcggaaccagaggcattgacaacagggttcgtagaagattcaagggccagatcttgatgcccaacattggttat
ggaagcaacaaaaaaacaaagcacatgctgccagtggttcgggaagtctctgtccacaacgtcaaggagctggaaagtgcgtgatgtgcaacaaatcttactgtgccgagatcgtcacaatgtttctccaagaaccgcaaagccatcgtggaagagctgcccaactggccatc
agagtcaccaacccaatgccaggctgcgcagtgaagaaatgagtag