

TAZ-tafazzin Gene

Reference: PTXBC011515

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BTHS; CMD3A; EFE; EFE2; G4.5; LVNCX; Taz1; tafazzin; protein G4.5		

For research use only.

GENE DESCRIPTION	tafazzin
GENE ID	6901
ORIGIN SPECIES	Human
ACCESSIONS	BC011515
VECTOR	pDONR223
NCBI SYMBOL	TAZ
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctctgcacgtgaagtggccgtcccccgggtgccgcctcacctggaccctggccagcagctcgtcatggccttggggacactacagctgcttctggaccagtgagtggtggccagggcggcagggcggcgggtaccatgcccgccggaggatcctgaaactc
cgccacatctggaacctgaagtgtgatgcttggacccctgcagctgcagacatctcttaccacagggagctacactccactcttcagcttgggcaagtgtgtgctgtgtgcccaggagatggcgtctaccagaagggatggacttcatttggagaagctcaacctggggactgggt
gcatacttcccagaaggaaatcggcgctgattgctgagtgcatctcaacccatcatcctgccccctgtgcatgtcgggtgagcctggggacggggacagagagatggcatctgggtggggcctgggactccctctggtcccaggctgccccctgctcccccacgtctgctctct
gtccactgtgtgcaggaaatgaatgacgtccttcaacagtcggccctacttcccccttggacagaaaatcactgtgtgatcgggaagccctcagtgccctgctgtactcagcggcctccggcggaacaagtggctgtggagatgcggaaagccctgacggacttcattca
agaggaattccagcatctgaagactcaggcagagcagctccacaaccacctccagcctgggagatag