

## SF3A3-splicing factor 3a, subunit 3, 60kDa Gene

Reference: PTXBC011523

|                |                                                                                                                                                                                                                      |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                  | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                 | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | PRP9; PRPF9; SAP61; SF3a60; splicing factor 3A subunit 3; SAP 61; pre-mRNA splicing factor SF3a (60kD); spliceosome associated protein 61; splicing factor 3a, subunit 3, 60kD; splicing factor 3a, subunit 3, 60kDa |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | splicing factor 3a, subunit 3, 60kDa                                                             |
| <b>GENE ID</b>          | 10946                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC011523                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | SF3A3                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggagacaatactggagcagcagcggcgtatcatgaggagaaggaacggctcatggacgtcatggctaaagagatgctcaccaagaagtcacgctccgggaccagatcaattctgatcaccgactcgggccaatgcaagataggtatgaggtcagtgaggaaacctgagggatt  
tgtatgatgataaggatgattacgaaaggaggagctcaatgccatttcaggacccaatgagtttgcgtgaattctataatagactcaagcaataaagggaattccaccggaagcaccacaatgagatctgtgtgccaatgctagtggaatttgaggaaacctctgaaggctcgagagaatccaa  
gtgaagaggcacaataacttggagggttcacagatgaggaggatattggtctgtatctcgtatccatgactgtttacctcaagtaacatgaaggcatctgagaagctggattatcacataacctgtccatcttgaccaattatttgacattctaaagaaaggaaatgcagagtataa  
gagatactagagatgctgttgcgtacgttcaggattacacagatagagtgaaagcctctccaagatcagaatgaacttttgggaagattcaggctgagtttgagaagaatgggagaatgggacatttctggatggccgaaagagacaagcagtgccctgacccatgctggagcccatc  
ttgacctctctgacattctctctcctggaggagttggcttctctgggtttggacagattgaaatctgctctcttagctttaggcttgaattgtggcgggaccctagaagagcagcccagagactatcagtagtaccagaaaggaaagtccttgagtcacttgataacctttgtttgccaataatcccaa  
gtcaaaaggccaccaagcgagacactgaaaggacaagacattgctttctagaagccagatctatgaatatgtagagattctcggggaacgcgacatctcactcatgaaatgtacagcgaagcaagccaggacagagagaagagcgaagaagaggaagaagagacatca  
gtgagagtgagagtgaaatgaagagaacgagatcatttacaaccccaaaacctgccatttggctgggatggcaaacctattccctactggtctgtataagcttcatggcctaataatcaactacaactgtgagatttgggaaactacacctaccgagggcccaagcctccagcgacact  
ttgctgaatggcgtcatgctcatggcatgaggtgtttgggcatcccaacactgctcatttctaatgtgacacagattgaagatgctctctctgtggcccaactgaaattgcagaaggttcagaacgatggcagcctgacactgaggaagaatatgaagactcaagtgggaattgtt  
gaataagaagacatacaggagatctgaaaagacaaggactgctctag