



Reference: PTXBC011568

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SIL1 nucleotide exchange factor; SIL1-like protein endoplasmic reticulum chaperone; SIL1 homolog, endoplasmic reticulum chaperone; nucleotide exchange factor SIL1; BAP; MSS; ULG5; BiP-associated protein		

For research use only.

GENE DESCRIPTION	SIL1 homolog, endoplasmic reticulum chaperone (<i>S. cerevisiae</i>)
GENE ID	64374
ORIGIN SPECIES	Human
ACCESSIONS	BC011568
VECTOR	pDONR223
NCBI SYMBOL	SIL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]