



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GTF3C6-general transcription factor IIIC, polypeptide 6, alpha 35kDa Gene

Reference: PTXBC011593

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf51; TFIIC35; bA397G5.3; general transcription factor 3C polypeptide 6; TFIIC 35 kDa subunit; general transcription factor IIC, polypeptide 6, alpha 35kDa; transcription factor IIC 35 kDa subunit; transcription factor IIC 35kDa; transcription fac		

For research use only.

GENE DESCRIPTION	general transcription factor IIIC, polypeptide 6, alpha 35kDa
GENE ID	112495
ORIGIN SPECIES	Human
ACCESSIONS	BC011593
VECTOR	pDONR223
NCBI SYMBOL	GTF3C6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcggcgcgagcagcgagtcaccagagcaggagaagacgaagaaaggaggagcagctgggtctgggtgaattatcaggaattatgattcagactctctcaaaatgtgaaaaaataatgcaagggtttgggcattgacactgagaggccattctgcaagtggacagctgtgctttctcgggagtatgaagacacctctagggacctgtgttatattggaagaaatgttgaacatgctgatacagaaggccaataataaaacagtgctaaataataaatgccatacaatgaagaagctcagcatgacagaagaactctctgacagagaagaagggaagggaagaaaaataggctgggtggaaatggctcgaataaagagataaagtctctcatcgcccaacatgatttgaactttctcatgataaataagaagcgaagaaggtgctgctctcttattcattacaccttcctgaagaaaggttctcttgaataaataacaaatgagctcccttagcgaatcgaatcgaatgaaatgaaattctctcttattcattacaccttcctgaagaaaggttctcttgaataaataacaaatgagctcccttag

More info: Please contact us for any question: contact@proteogenix.fr