



Reference: PTXBC011625

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NOLZ-1; NOLZ1; Nlz2; zinc finger protein 503		
For research use only.			

GENE DESCRIPTION	zinc finger protein 503
GENE ID	84858
ORIGIN SPECIES	Human
ACCESSIONS	BC011625
VECTOR	pDONR223
NCBI SYMBOL	ZNF503
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr