



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

BMI1-BMI1 polycomb ring finger oncogene Gene

Reference: PTXBC011652

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BMI1 proto-oncogene, polycomb ring finger; polycomb group protein Bmi1; BMI1 polycomb ring finger proto-oncogene; BMI1 polycomb ring finger oncogene; FLVI2/BMI1; PCGF4; RNF51; flvi-2/bmi-1; polycomb complex protein BMI-1; B lymphoma Mo-MLV insertion regio		

For research use only.

GENE DESCRIPTION	BMI1 polycomb ring finger oncogene
GENE ID	648
ORIGIN SPECIES	Human
ACCESSIONS	BC011652
VECTOR	pDONR223
NCBI SYMBOL	BMI1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcatcgaacaacgagaatcaagatcactgagctaaatccccactgtagtggctgcttggggagggtacttcattgatgccacaaccataatagaatgtctactctctctgtaaaacgtgtattgtctgtacctggagaccgaagtatgtctatttggatgtccaagtgcacaagacc
actgactgactgaataataagcttcagataaaactcccaagatgatatacaaaattgtccaggctgtgtcccaaaatgaatgaagagaagaagggtattgtcagctcatctctctgtgatctgtcccaatgctccttaagatagaatagatgagatgagaataagagaataata
acgtgatgatgataataagcttctcattcgaattcttcaggacagaatgatgagtcagaagaataaaccaagaagaacgaatctaaaggagtggaatgaataaagaatacttaccatggtcccgacgaacgaatgatgtgatgacttaagaagatttcagatgagaataaagtacatcata
atacttccagattgatgcatgatgaggagggaacctttaaaggattattatacctaattgatatctctacatttatactctggagaaggaaatggtccacttccattgaataacagatgtgacctactgtaaaagaatgaagatcagtcaccagagagatggactgacaatgctggagaact
ggaagtgaaatctggaggatgacaaggaccccgacaggagattcttccctcaactctctctgtttgctagccccgactccagtgcatgtctctcatccacagtttctccatatttccagtagtatgaatggaaccagcaacagccccgcggtgaaccaccaatctcttttgcgaata
gacctcgaaactcatcgaataaagcttcacgacaaactctgttgta

More info: Please contact us for any question: contact@proteogenix.fr