



Reference: PTXBC011669

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DCE; Nbla03646; SELADIN1; seladin-1; delta(24)-sterol reductase; 3 beta-hydroxysterol delta 24-reductase; desmosterol-to-cholesterol enzyme; diminuto/dwarf1 homolog; seladin 1; selective AD indicator 1; 24-dehydrocholesterol reductase		

For research use only.

GENE DESCRIPTION	24-dehydrocholesterol reductase
GENE ID	1718
ORIGIN SPECIES	Human
ACCESSIONS	BC011669
VECTOR	pDONR223
NCBI SYMBOL	DHCR24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]