

TTC19-tetratricopeptide repeat domain 19 Gene

Reference: PTXBC011698

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	2010204O13Rik; MC3DN2; tetratricopeptide repeat protein 19, mitochondrial; TPR repeat protein 19; tetratricopeptide repeat domain 19		

For research use only.

GENE DESCRIPTION	tetratricopeptide repeat domain 19
GENE ID	54902
ORIGIN SPECIES	Human
ACCESSIONS	BC011698
VECTOR	pDONR223
NCBI SYMBOL	TTC19
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagatgagccagaagaggtgagtttaatttgcagcgtctctgctcgcctatcagactgataacaagaaggccatcacttacattgatgttgcccaacttagcatttatacggggtcagcttgaataatgctgaacaacttttaaagcaaatgagttacctctggagggg
gcatgaagcaggaggacaatgcaataattgaaattccctaaagctggccatgatctatgctgcgcagacagacaggaatttgcgttgctgctgaattctgcatttcaactctagaggaataatgaaagagaaaggaattagcagaagacattatgtcagtggaagagaaagcca
ataccacctctcttggcagtgcttagacgcctgtgctcgtacctctgttctccaagcagccgtcacaggcacaaaggatgtatgaaaagctctgcagatttctgaagaatacaaggagaagacaccacagaccattgtgctgatgagtgacctggctactacctggatgcac
agggccctttgatgagccctatttatgcaaaaggcatcagatctggcaagacagataaatcatcctgagctacacatggtactcagtaactagctgcagttttgatgcacagagaacgatatacacaagcaaaagagatctaccaggaagcactgaagcaagcaagctgaaaaaa
gatgaaatttctgatacacatcagggaagagttggctgagctgtcaagaaaagtagaccttggacaaattctgtcaagctctaa