



Reference: PTXBC011726

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GEF1; LBCL2; LSC; P115-RHOGEF; SUB1.5; rho guanine nucleotide exchange factor 1; 115 kDa guanine nucleotide exchange factor; 115-kD protein; Lsc homolog; Rho guanine nucleotide exchange factor (GEF) 1; p115RhoGEF		

For research use only.

GENE DESCRIPTION	Rho guanine nucleotide exchange factor (GEF) 1
GENE ID	9138
ORIGIN SPECIES	Human
ACCESSIONS	BC011726
VECTOR	pDONR223
NCBI SYMBOL	ARHGEF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1