

RPP21-ribonuclease P/MRP 21kDa subunit Gene

Reference: PTXBC011730

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf135; CAT60; ribonuclease P protein subunit p21; RNaseP protein p21; ribonuclease P/MRP 21kDa subunit; ribonucleoprotein V; ribonuclease P/MRP subunit p21		

For research use only.

GENE DESCRIPTION	ribonuclease P/MRP 21kDa subunit
GENE ID	79897
ORIGIN SPECIES	Human
ACCESSIONS	BC011730
VECTOR	pDONR223
NCBI SYMBOL	RPP21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggggccgggtgaaggaccgcgagggccttcagaggctcaacttcctgtaccaggccgccattgtctctgccagaccgcgagaccgagaaccaggcgtggcgaggtttactgtctacactgagaggaccattgcgaagcggctcgtcttgcggcgggacccctcgggtgaagagga
ctctctgtcggaggtcctcttccctcctcgtccgggctcactgcacccaccgacagagacgtgcaggggacagcgtggaccgtacagactgcctaacatgccagcgcagccaacgttctcctaatgatccgggcatttactctggggagacaggcctgagggccagctcggg
agccaagcagattccaaccactacaaccttgccaacacagcccactccatttcagaccgcttctgaggagaaaatgcagactcagggttcagtaaccagtga