

## IMP3-IMP3, U3 small nucleolar ribonucleoprotein, homolog (yeast) Gene

Reference: PTXBC011745

|                |                                                                                                                                                                                                                                 |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | IMP3, U3 small nucleolar ribonucleoprotein; U3 snoRNP protein IMP3; IMP3, U3 small nucleolar ribonucleoprotein, homolog; U3 small nucleolar ribonucleoprotein protein IMP3; BRMS2; C15orf12; MRPS4; U3 snoRNP protein 3 homolog |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | IMP3, U3 small nucleolar ribonucleoprotein, homolog (yeast)                                       |
| <b>GENE ID</b>          | 55272                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC011745                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | IMP3                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggtgcggaagcttaagttccacgagcagaagctgctgaagcaggtggacttcgaactggaggtcaccgaccacaacctgcacgagctgcgctgctcggttaccggctgcagcggggaggactacacgcgctacaaccagctgagccgtgccgtgcgtgagctggcgcggcgctgcgcgactgcccgaacgcgaccagttccgctgcgcgcttcggcgctctggacaagctgtatgctctcggcttgggtgccacgcgcggttcgctggagctctgcgacttcgtcagcgctcgtcctctcgcggccgctccccaccgtgctcctcaagctgcgatggcgagcagcacttcaggctgcccgtggccttctggagcaaggacgctacgctgggcccctgacgtggttaccgaccccgcttctgtcagcgcgagcatggaggactttgcaacttgggtggactcgtccaagatcaagcggcacgtgctagagtacaafgaggagcgcgatgacttcgatctggaagcctag