

C16orf42-chromosome 16 open reading frame 42 Gene

Reference: PTXBC011785

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16orf42; ribosome biogenesis protein TSR3 homolog; TSR3, 20S rRNA accumulation, homolog; TSR3, acp transferase ribosome maturation factor		

For research use only.

GENE DESCRIPTION	chromosome 16 open reading frame 42
GENE ID	115939
ORIGIN SPECIES	Human
ACCESSIONS	BC011785
VECTOR	pDONR223
NCBI SYMBOL	C16orf42
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgaggagccacttgcctgttgcctacctgggtggccgcaaccccgtagaactatggccgaccttctgcgtggaagcgtttgctgccaccttctgcatcgtaggcttccagaccttgcgtcatctttgtgcggaagttaaatggggcaagggtcttggacctgaa
ccgccagctcctggacaagtacggcctgcggcagccggaggaggtgctgcaggcggagcaggagttcttggccaatgccaaggagagccccaggaggaggagatcgtaccccttcgatgtggatcaggagagagatttggaaacccaacaggcctgtggccagacccg
gctgccctcggacactgatgacagtgtgctgtgaggaccaggcctggcggcggagcggaggagccagcagcgtgctgtgaagaggagcagacgcaggacggggggctgaggccaggcccgctgagggttggaaaggaaatcaagaacggcagagagactga