

DNAJB12-DnaJ (Hsp40) homolog, subfamily B, member 12 Gene

Reference: PTXBC011812

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DJ10; dnaJ homolog subfamily B member 12; DnaJ (Hsp40) homolog, subfamily B, member 12; DnaJ heat shock protein family (Hsp40) member B12		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily B, member 12
GENE ID	54788
ORIGIN SPECIES	Human
ACCESSIONS	BC011812
VECTOR	pDONR223
NCBI SYMBOL	DNAJB12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctttggcggcggttccttctagtaacgtccacgtctacagcaacggccgcatgcgtatactaccagcaaaaggcaggaccgcagggacaaccagggtgatggcgggtaggggtgtttgtgcagctgatgcctatcctatcctgattctcgtgcagctctcagccagctcatg
gtctccagtcaccctacagctctgagtcgaagaccgtccgtggccacatccacagcgagtcactgaccacctgggtgtctctactatgtgggagacacttttccgaagagtacacaggctccagcctcaaaacagtcgagcgggaatgtggaagatgattatcgccaaacctccgga
acaactgtggaaggagaagcagcagaaggaggcttgcgtgaccgggcagcgtactttggcgacacagatatgtaccacagagcacagaagatgggcacccccagctgcagccgactgtcagaggtgcaggcctccctgcatggatag