



Antibodies and Proteins

## ProteoGenix SAS

## Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

[www.ProteoGenix.science](http://www.ProteoGenix.science)

## SLC48A1-solute carrier family 48 (heme transporter), member 1 Gene

Reference: PTXBC011823

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

### Synonyms

HRG-1; HRG1; hHRG-1; heme transporter HRG1; heme-responsive gene 1 protein homolog; solute carrier family 48 (heme transporter), member 1; solute carrier family 48 member 1

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | solute carrier family 48 (heme transporter), member 1                                             |
| <b>GENE ID</b>          | 55652                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC011823                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | SLC48A1                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

[illegible]