

HEL308-DNA helicase HEL308 Gene

Reference: PTXBC011863

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DNA helicase HEL308; HEL308; helicase POLQ-like; POLQ-like helicase; mus308-like helicase; helicase, POLQ-like		

For research use only.

GENE DESCRIPTION	DNA helicase HEL308
GENE ID	113510
ORIGIN SPECIES	Human
ACCESSIONS	BC011863
VECTOR	pDONR223
NCBI SYMBOL	HEL308
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgaatgtggtcccgcatccgcccgggtgtctctcccaaaaggaaccgtccaagcttgggtgtattttggcgctccaccgcccggcgagctcgagccggagatgaggggaaagaggaggaggaatgtggtgctgagaacaggaggcggaaaccgcggcggtact
gcccgttgagggtacagcccccttctactctcagattccccggaatgtctcgtcttggagggtggtgatacaaacccgacctcctacgtcacatgccactgacagagggtgggagaccagcctaataatgacagtgaagtggacatgtttggtgactatgatacttactgaaactccttata
gctcaagtgtgacgacctggaacaaaaatataatgcaactccctgaacataagaacatgctacagactttgccactgaaaactcttctcggaagtatacaaaaacacacagcattactaccataggcaaccttactgaattacaactgataagcacacagagaaccagagtgatagaag
gtgtcactattgaacctggagctgactttgtatgatgtacctctcacaggctatatcttgaataatgacagaactcttcaaatgatttgggtgatcttctatgaaagaaaggattggaagtcctctcacaacactgtgaatgaggaaactgcccataattgcatagagcaacccagc
aaaatgatgagctcttccaaagtcagaactagttcagatatgaacaggagaaaaagtattaaagatcatctaaaaatgccatga