

FKBP9L-FK506 binding protein 9-like Gene

Reference: PTXBC011872

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HMG20; polyubiquitin-C; ubiquitin C		

For research use only.

GENE DESCRIPTION	FK506 binding protein 9-like
GENE ID	360132
ORIGIN SPECIES	Human
ACCESSIONS	BC011872
VECTOR	pDONR223
NCBI SYMBOL	FKBP9L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacatgggtctcagagagatgtgcgttggcgagaaacggacagtatcattccgctcacctgggctatggggaagctggcgtggatggagaagtggcggcagtgccgtattagtgttgacattgagctgctggagctgggtggcctcccgaggggtacatgttcatatgga
atggtaggtgtcacccaaccttttgaagaaattgacaaggatggcaacggagaaagtcctcctggaagagttctcagagtacattcacgccaggtggcatctggcaagggaactcgtcctggctttgatgctgagctgattgtgaagaatatgttccaaccaggaccggaatgga
gatgggaaggtcacagctgaggaatttaactcaagaccaggaagccaacaggtgaactctaa