

C1orf212-chromosome 1 open reading frame 212 Gene

Reference: PTXBC011880

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	alternative protein C1orf212; UPF0767 protein C1orf212; C1orf212; small integral membrane protein 12		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 212
GENE ID	113444
ORIGIN SPECIES	Human
ACCESSIONS	BC011880
VECTOR	pDONR223
NCBI SYMBOL	C1orf212
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggcctgtgttttgaccgtgggttcgtacatgctcctatgtcacattccctgttgcctcgtggtcgggctgtgggttaccacctgggaatggttcacagggaagacccccagcccgaggaggaaaagagcatctcagagcgggggaggatcgcaagctggatgagctt
ctaggcaaggaccacacgcaggtggtgagccttaaggacaagctagaattgccccgaaagctgtgctgaacagaaacgccagagaagaattaa