

COQ4-coenzyme Q4 homolog (*S. cerevisiae*) Gene

Reference: PTXBC011895

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ubiquinone biosynthesis protein COQ4 homolog, mitochondrial; CGI-92; COQ10D7; coenzyme Q biosynthesis protein 4 homolog; coenzyme Q4 homolog; coenzyme Q4		

For research use only.

GENE DESCRIPTION	coenzyme Q4 homolog (<i>S. cerevisiae</i>)
GENE ID	51117
ORIGIN SPECIES	Human
ACCESSIONS	BC011895
VECTOR	pDONR223
NCBI SYMBOL	COQ4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgactctgctgcgccctgtctcctcgtcggctctcggggtcccggtacagcggcctgcggcagaaatgccctccgggctagagcagcggcgccggcccgctatactcgaccacctccccacctccccgctcagaaagcgtgttggccgccgctccgcggcgat
ggcgctctataacccctaccgccagacatggctcagcttctaggggagaccacagcagcaccgacccctgaaaggctcctcagggaccagatgaggaggatccagagggtgccagatcctgcaggagcgtccccggaattcgacatccacctcgacctgggcaagctccagagcct
gccggaaggctccctcggctcgcgagatctcctgttctggtgagtgacaggggtctccccagacaccgagcaccaccgctcgtggatgatgaggagctagcgtatgtattcagcgggtaccgggaggtgcacgacatgcttcacacctgctgggagatgccaccaacatcctggg
ggagatcgtggtgaaatggttgaggctgtccagactggcctgccatgtgcctcctgggtgcattctttggaccgacatccgacttggcgtcagagcctgcaagtgcgtcgtcggagttgatccatgggcccgttcagaacggcgagagccccatgtgtcctcaacctgtactatgagc
ggcgctgggagcagtcctgagggtcgtcgggagggagctgggcattacagcaccacccatgcacgtccagggttggcctga