

MYH9-myosin, heavy chain 9, non-muscle Gene

Reference: PTXBC011915

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BDPLT6; DFNA17; EPSTS; FTNS; MHA; NMHC-II-A; NMMHC-IIA; NMMHCA; myosin-9; cellular myosin heavy chain, type A; myosin, heavy chain 9, non-muscle; non-muscle myosin heavy chain 9; non-muscle myosin heavy chain A; non-muscle myosin heavy chain IIa; non-musc		

For research use only.

GENE DESCRIPTION	myosin, heavy chain 9, non-muscle
GENE ID	4627
ORIGIN SPECIES	Human
ACCESSIONS	BC011915
VECTOR	pDONR223
NCBI SYMBOL	MYH9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaaacggatgggctgagggcaaacctgccgaataagcctctctcctgcagcctgagatggatggagacagacacacagcctccccctccagaccccgagcagcctctccccacctcttgggactgctgtgaacatgcctctcctgcctccgccccgtccccca
tcccgtttccctccagggtgtgttgagggcatttggcttctctgctgcacccctccagctccctccctgctcagaatctgatacacaagagacagggcccgggccagcgagagcgaccagcaggtctctcagcctctcttgcacaaaagcacaagatgttgaggcgagcagg
caggccccggggagggccagagttttctatgaatctattttcttcagactgaggccttttgtag