



Reference: PTXBC011948

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PL6; transmembrane protein 115; PP6; placental protein 6; testicular tissue protein Li 201		

For research use only.

GENE DESCRIPTION	transmembrane protein 115
GENE ID	11070
ORIGIN SPECIES	Human
ACCESSIONS	BC011948
VECTOR	pDONR223
NCBI SYMBOL	TMEM115
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1