

WDR32-WD repeat domain 32 Gene

Reference: PTXBC011959

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WDR32; DDB1- and CUL4-associated factor 10; WD repeat domain 32; WD repeat-containing protein 32; DDB1 and CUL4 associated factor 10		

For research use only.

GENE DESCRIPTION	WD repeat domain 32
GENE ID	79269
ORIGIN SPECIES	Human
ACCESSIONS	BC011959
VECTOR	pDONR223
NCBI SYMBOL	WDR32
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgaatgagggttaaccaccagattgtccaaatgttgattcaacgtcctctggatctcttaattttgcatgaccttgacttaactaagctttagaagtaggcagctatccatttaagagcaagaagaactactcaagtcagatttgaccacttcatcaagttcatctggtcctagagttctg
gctcaccttgatcatagtgattctaatctctgagaacacatgacagcctctcaagagaaggagtttcaccacgaaatagcttggaagtcgtaacccctgaagttcttggtgagagtgaccacgaaactgcatcagctccttacagctgcacccaaaggctgggccactctcttc
gggtctcaagcaacagtgatgatgaggagtgacttggtctatgaattccaagaaggagctccagtgcgacctgtctcacctcgtgttctctacgactgactcattacattgaagaagccaatgtaggcaggggttacatcaagaacttgcttcagccccgatggacgaatgatttctcc
cacatggctatgggattcgttgggatttgacaacagtgacgtgaactgttgactgctgccccaaagaagccagtcctcgcgggtgatccgttctgtactctcataatgatgtgtactgacacaaagtctctccaacacattgtcagattgctcaggtgcttagtgacgggt
ttctttgatcagccaaagttag