

PSMD10-proteasome (prosome, macropain) 26S subunit, non-ATPase, 10 Gene

Reference: PTXBC011960

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dJ889N15.2; p28; p28(GANK); 26S proteasome non-ATPase regulatory subunit 10; 26S proteasome regulatory subunit p28; ankyrin repeat protein; hepatocellular carcinoma-associated protein p28-II; proteasome (prosome, macropain) 26S subunit, non-ATPase, 10; pr		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) 26S subunit, non-ATPase, 10
GENE ID	5716
ORIGIN SPECIES	Human
ACCESSIONS	BC011960
VECTOR	pDONR223
NCBI SYMBOL	PSMD10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggggtgtgtctaacctaattggtctgcaacctggcctacagcgggaagctggaagagtggaaggagattctgcccataaatccctggctactagaactgaccaggacagcagaactgcattgcactgggcatgctcagctggacatacagaattgtgaattttgtgcaa
cttggagtgccagtgaatgataaagacgatgcaggttggtctcctctcatattgcggctctgctgcccgggatgagattgtaaaagccctctgggaaaaggctcgaagtgaatgctgcaatcaaatggctgtactccttacattatgcagcttcgaaaaacaggcatgagatcgtgtc
atgttactggaaggcgggctaataccagatgctaaggaccattatgagctacagcaatgcacgggcagcagccaagggttaacttgaagatgattcatatccttctgtactacaaagcatccacaacatccaagacactgagggttaacactcctacacttagcctgtgaggagag
agtggaaagcaaaactgctggtgtcccaaggagcaagtaattacattgagaataaagaagaaagacaccctgcaagtggccaaaggctggcctgggttaataactcaagagaatggtggaaggttaa